



B"H
Maxima Medical LTD
Haogen 10, Northern Industrial park
POB 4124, Ashdod 77141

16th of Iyar 5772, 08/05/12

Shabbat Certification

We hereby certify

that the CPAP devices sold by **Maxima Medical**
and manufactured by **BMC models**
Resmart CPAP / Resmart Auto CPAP
For the treatment of sleep apnea
can be used on the Jewish Sabbath and Holidays
according to the following stipulations:

A. Programming (before Shabbat):

1. Program the display to be at a fixed lighting (On or Off) and not automatic mode (Auto Light).
2. **We recommend that the machine be connected to a Shabbat clock** which was set before the beginning of Shabbat. To do so, it's necessary to program the functions **Auto-off and Auto-on mode to disable**.
3. If you use a Shabbat clock, the **buzzer** should be disconnected (by company's technicians) so that it won't beep when the power is turned off by the clock, till it finishes the internal battery.
4. If you use a Shabbat clock, it is recommended to program the mode **Alert Disable** to prevent beeping when taking off the mask after sleeping on Shabbat.

B. Using Shabbat :

1. Do not press any button on Shabbat (except in an emergency and with '*Shinuy*' only).
2. It is recommended to cover the display and not to pay attention to the information displayed there in.
3. As mentioned, **it is recommended** to use with a Shabbat clock, and make the guidelines above Section A 2,3,4 .
4. If necessary (due to breathing difficulties, severe discomfort, etc.) **it's allowed to operate the device by breathing into it and automatically disconnect in the absence of breath**. This method does not require a Shabbat clock, instead turn the device on before Shabbat.
In this situation, program the functions **Auto off and Auto-on to mode Enable**.
With this method too, you should program **Alert Disable** .
5. On Shabbat you can use the moisture facility. It's advisable to fill the container before Shabbat.

With Blessings for a restful and healthy Shabbat,

Rabbi Yisrael Rozen, eng.

Dean, The Zomet Institute

\\cpap-maxima-eng-f.docx\רפואי ושמייעד\1\wordocs\israel\technology\yishurim\